

Continuous Passive Motion is a term often used in rehabilitation. The idea is that if you keep a joint moving you can decrease stiffness and increase range of motion. Golden's AutoFlex power lumbar system was created with that idea in mind.

- AutoFlex is a metal, scissor-like mechanism, not your traditional lumbar air bladder. This ensures durability and does not run the risk of deflating.
- The continuous passive motion can provide relief to a sore back by stretching the spine and lower muscles
- AutoFlex can be used as stationary lumbar support that can be adjusted to your comfort level.
- Three cycles with three different intensities are provided and can be used based on need

Autostart

Autostop

AutoFlex

• Best if used in Zero-gravity or Trendelenburg* position

*Consult doctor before experiencing the Trendelenburg position.

Decrease Stationary Lumbar Support

Increases Stationary Lumbar Support

Give relief to the lumbar with



Available on select models.



SHT-UCAFLEX.010318